NEW LIFE SPECTRUM®

High-Density Color Enhancing Gel Mix


Mixing Instructions





After introducing Gel, observe your aquarium. Food should be substantially consumed within approx. 15 min. Adjust amount of Gel mixed as needed

Remove uneaten portion within 4-6 hours.



Preparation Suggestion:

You can make a larger portion than your animals can eat in one feeding and divide (cut) the Gel into multiple servings to save for later.




Prepared Gel Storage:

Freezer: Will preserve indefinitely.
Seal in freezer bag to prevent freezer burn.

Refrigerator: Keeps 7 days. Store in container covered or sealed with plastic wrap.

Feeding/Training Tip:

Can spread or stick gel onto surfaces in the aquarium, such as decorative rocks, coral skeletons etc. — useful for training those used to foraging and feeding on rocks and corals.







ADDITIONAL TIPS

Simple Food Cubes:

If you want small cube portions great for many uses and similar to frozen food cubes, you can pour the heated mix into an ice tray. Party ice trays are perfect, as are silicon ice trays that allow you to "flex" the tray and pop-out the cubes.

Feed Training Station:

Because Nutri/Gel is a mix, there's an even better way to make a "training station" in your tank than spreading or sticking the gel onto a rock or coral. Instead, pour the hot mix onto the item before it cools, allowing it to coat the item. Do not put too much, and start with a conservative amount.

NOTES

 All aquarium communities are unique. When using Nutri/Gel, adjust feeding quantities based on your observation of the tank.

 This product should allow for multiple preparations/ servings per container & comes in multiple jar sizes.



31950 SW 197 Ave., Homestead Florida (USA) 33030 www.nlsfishfood.com